


## smart moves complete with ComfortTrack™ material

Made exclusively for clear aligners, this revolutionary material's elastomeric properties generate an ideal initial force and maintain it over time. **The result? Superior patient comfort and more effective aligners.**

### ComfortTrack™ Material vs. Other Clear Aligner Materials



\*Force Retention at 37 Degrees Celsius

Compared to standard aligner plastics, ComfortTrack's multi-layered structure:

- » Generates 40% less initial force for reduced patient pain
- » Retains 200+% greater force on day 7

**You'll feel the difference,  
and your patients will too.**

Get in touch with us to learn more from a smart moves expert.

800.828.7626  
smartmovesaligners.com



Great Lakes™  
DENTAL TECHNOLOGIES